

1) Measure the DNA concentration (ng/μl) of each assembly piece.

2) Add 100 ng of the linearized vector backbone and equimolar amounts of the other assembly pieces to a 15 μl total volume assembly reaction mixture as follows:

	linearized vector backbone (100 ng)
+	each additional assembly piece (to equimolar with backbone)
+ 1.5	μl 10X NEB T4 Buffer
+ 0.15	μl 100X BSA*
+ 1	μl <i>Bsal</i>
+ 1	μl NEB T4 Ligase, 2 million cohesive end units / mL
+ _____	dH ₂ O to
15	μl

NOTE: It is **essential** to use a High Concentration Ligase

* *Bsal* is only 10% active at 37 C without the addition of BSA.

3) Perform the assembly reaction in a thermocycler as follows:

either (following Engler 2009):

3	min @ 37 C }
4	min @ 16 C } 25 cycles

5	min @ 50 C }
5	min @ 80 C } 1 cycle

or, alternatively (modified from Engler 2008)

1	hour @ 37 C	1 cycle
---	-------------	---------

5	min @ 50 C }
5	min @ 80 C } 1 cycle

NOTE: If any of the assembly pieces contain an internal *Bsal* site(s), it would be first preferable to silence the internal *Bsal* site(s) through point mutation(s); a second option, if the internal *Bsal* site overhang(s) are not cohesive with the other assembly overhangs, is to adjust the thermocycling parameters to terminate after a ligation step (e.g. skip the final cycle at 50 and then 80 C).

4) Transform 5 μl of the assembly reaction into 100 μl of competent *E. coli* and/or run a diagnostic agarose gel to check for successful assembly.