

1) Measure the DNA concentration (ng/μl) of each assembly piece.

2) Add 100 ng of the linearized vector backbone and equimolar amounts of the other assembly pieces to a 25 μl total volume assembly reaction mixture as follows:

		linearized vector backbone (100 ng)
+		each additional assembly piece (to equimolar with backbone)
+ 5	μl	5X HF Phusion Reaction Buffer
+ 1	μl	10 mM dNTPs
+ 0.75	μl	DMSO
+ 0.5	μl	2U/μl Phusion Polymerase
+ _____		dH ₂ O to
25	μl	

3) Perform the assembly reaction in a thermocycler as follows:

	30	sec	@ 98 C	1	cycle
	10	sec	@ 98 C }		
	30	sec	@ 55 C }	1 to 15 cycle(s)**	
length* (kb) x 15	sec	@ 72 C }			
	10	min	@ 72 C	1	cycle

*The total length of the assembled product (in kb)

**The number of repeated cycles should exceed the number of assembly pieces

4) Transform 5 μl of the assembly reaction into 100 μl of competent *E. coli* and/or run a diagnostic agarose gel to check for successful assembly.