

SLiCE Reaction

1) Prepare 10X SLiCE Buffer in a 1.5 mL tube :

500 μ L 1 M Tris-HCl pH 7.5
+ 50 μ L 2 M $MgCl_2$
+ 100 μ L 100 mM ATP
+ 10 μ L 1 M DTT
+ _____ ddH₂O to
1 mL

Store at -20 C in 40-60 μ l aliquots.

2) Add the following ingredients into a 0.2 mL tube in this order and **vortex**:

linearized vector backbone (50-200 ng)
+ each additional assembly piece (1:1-1:10 molar ratio of
vector:insert)
+ 1 μ L 10X SLiCE buffer
+ 1 μ L PPY SLiCE extract
+ _____ ddH₂O to
10 μ L

3) Incubate the SLiCE reaction mix as above at 37 C for 1 hour using a PCR machine or water bath, and then place on ice.

4) Transform 1 μ L of the assembly reaction into 100 μ L of competent *E. coli* and/or run a diagnostic agarose gel to check for successful assembly.

5) For electroporation, transform 1 μ L into 20 μ L electrocompetent cells. For large recombinant DNA, electroporation is required. In complex cloning, electroporation is recommended, as it is 10-100 times as efficient as chemical transformation.